

IL13

Research Grade Lebrikizumab

Catalog No.	DHE07702A	Quantity:	100 µg
	DHE07702B		1.0 mg

Alternate Names: Interleukin-13, IL-13, NC30, MILR1444A, PRO301444

Description: Lebrikizumab is a humanized monoclonal antibody that targets IL-13. It is under development for the treatment of asthma and atopic dermatitis.

Concentration: 1.0 mg/ml, lot specific

UniProt ID (target): P35225

Source: XtenCHO cell line

Isotype: Humanized IgG4 kappa

INN: Lebrikizumab

Specificity: Human IL-13

Endotoxin: ≤ 0.01 EU/µg by LAL analysis

Formulation: 0.01M PBS, pH 7.4

Purification: Protein A/G affinity chromatography

Purity: > 95% as determined by SDS-PAGE

Storage & Stability: Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. **Avoid freeze/thaw cycles.**

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.