

IL13, IL4

Research Grade Romilkimab

Catalog No.	DHE07704A DHE07704B	Quantity:	100 µg 1.0 mg
Alternate Names:	Interleukin-13, IL-13, Interleukin-4, IL-4, SAR156597, huTBTI3_2_1		
Description:	Romilkimab is a chimeric/humanized bispecific antibody that targets IL-13 and IL-4		
Concentration:	1.0 mg/ml, lot specific		
UniProt ID (target):	P35225, P05112		
Source:	XtenCHO cell line		
Isotype:	Bispecific IgG4, [VH-VH'-H-Gamma1-VL-VL'-C-kappa]-dimer		
INN:	Romilkimab		
Specificity:	Human IL-13 and Human IL-4		
Endotoxin:	≤ 0.01 EU/µg by LAL analysis		
Formulation:	0.01M PBS, pH 7.4		
Purification:	Protein A/G affinity chromatography		
Purity:	> 95% as determined by SDS-PAGE		
Storage & Stability:	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.