

IL13

Research Grade Abrezekimab

Catalog No.	DHE07705A	Quantity:	100 µg
	DHE07705B		1.0 mg

Alternate Names: Interleukin-13, IL-13, UCB-4144, CDP7766, VR-942

Description: Abrezekimab is a humanized monoclonal antibody that targets IL-13 in development for asthma. It was designed as a high-affinity antibody fragment suitable for inhalation, to provide a more rapid onset of action at lower doses than subcutaneous systemic administration.

Concentration: 1.0 mg/ml, lot specific

UniProt ID (target): P35225

Source: XtenCHO cell line

Isotype: Human Fab-G1 kappa

INN: Abrezekimab

Specificity: Human IL-13

Endotoxin: ≤ 0.01 EU/µg by LAL analysis

Formulation: 0.01M PBS, pH 7.4

Purification: Protein A/G affinity chromatography

Purity: > 95% as determined by SDS-PAGE

Storage & Stability: Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. **Avoid freeze/thaw cycles.**

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

