

IL13

Research Grade Dectrekumab

Catalog No.	DHE07707A DHE07707B	Quantity:	100 µg 1.0 mg
Alternate Names:	Interleukin-13, IL-13, QAX-576		
Description:	Dectrekumab is a fully human investigational anti-interleukin 13 monoclonal antibody, with anti-inflammatory potential. It was developed for asthma and seasonal allergic rhinitis.		
Concentration:	1.0 mg/ml, lot specific		
UniProt ID (target):	P35225		
Source:	XtenCHO cell line		
Isotype:	Human IgG1 kappa		
INN:	Dectrkumab		
Specificity:	Human IL-13		
Endotoxin:	≤ 0.01 EU/µg by LAL analysis		
Formulation:	0.01M PBS, pH 7.4		
Purification:	Protein A/G affinity chromatography		
Purity:	> 95% as determined by SDS-PAGE		
Storage & Stability:	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

