

## KDR

### Research Grade Tanibirumab

<b>Catalog No.</b>	DHE16604A DHE16604B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	Olinvacimab, SSS-23, TTAC-0001, tanibirumab, CAS: 2095504-49-5		
<b>Description:</b>	Research Grade Tanibirumab is a human monoclonal antibody which targets human Vascular endothelial growth factor receptor 2 (CD309, FLK-1)		
<b>UniProt ID (target):</b>	P35968		
<b>Target:</b>	Human VEGFR-2 (CD309, FLK-1)		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Source:</b>	XtenCHO		
<b>Isotype:</b>	Human IgG1 kappa-lambda		
<b>Formulation:</b>	0.01M PBS buffer, pH 7.4		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Endotoxin Level:</b>	≤ 0.01 EU/µg by LAL analysis		
<b>Purification:</b>	Protein A/G affinity chromatography		
<b>INN:</b>	Tanibirumab		
<b>Applications:</b>	Functional studies		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. <b>Avoid freeze/thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



**Cell Sciences**<sup>®</sup>  
65 Parker Street  
Unit 11  
Newburyport, MA 01950

Toll Free: 888-769-1246  
Phone: 978-572-1070  
Fax: 978-992-0298

E-mail: [info@cellsciences.com](mailto:info@cellsciences.com)  
Website: [www.cellsciences.com](http://www.cellsciences.com)