

SNCA

Research Grade Anti-Human SNCA/Alpha-synuclein (Lu AF82422)

Catalog No.	DHE22008A DHE22008B	Quantity:	100 µg 1.0 mg
Alternate Names:	Alpha-synuclein, Non-A beta component of AD amyloid, Non-A4 component of amyloid precursor, NACP, PARK1, SNCA, Lu AF82422		
Description:	Research Grade Anti-Human SNCA/Alpha-synuclein (Lu AF82422) is a human monoclonal antibody (mAb) that recognizes all major species of alpha-synuclein (monomeric and oligomeric; N- or C-terminal truncated forms).		
UniProt ID (target):	P37840		
Specificity:	Human Alpha-synuclein		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO cell line		
Isotype:	Human IgG1		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com