

CD200

Research Grade Samalizumab

Catalog No.	DHE30301A DHE30301B	Quantity:	100 µg 1.0 mg
Alternate Names:	ALXN6000, Anti-CD200, hB7V3V2 G2G4, CAS: 1073059-33-2		
Description:	Research Grade Samalizumab is a humanized monoclonal antibody which targets human OX-2 membrane glycoprotein		
UniProt ID (target):	P41217		
Specificity:	Human OX-2 membrane glycoprotein (CD200)		
Source:	XtenCHO cell line		
Isotype:	Human IgG2-IgG4 kappa		
INN:	Samalizumab		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	0.01M PBS, pH 7.4		
Purity:	> 95% as determined by SDS-PAGE		
Purification:	Protein A/G affinity chromatography		
Endotoxin:	≤ 0.01 EU/µg by LAL analysis		
Storage & Stability:	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

