

CCR2

Research Grade Plozalizumab

Catalog No.	DHE33201A DHE33201B	Quantity:	100 µg 1.0 mg
Alternate Names:	MLN-1202, anti-CCR2, hu1D9, CAS: 1610761-46-0		
Description:	Research Grade Plozalizumab is a humanized monoclonal antibody which targets human C-C chemokine receptor type 2		
UniProt ID (target):	P41597		
Specificity:	Human C-C chemokine receptor type 2 (CD192)		
Source:	XtenCHO cell line		
Isotype:	Human IgG1 kappa		
INN:	Plozalizumab		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	0.01M PBS, pH 7.4		
Purity:	> 95% as determined by SDS-PAGE		
Purification:	Protein A/G affinity chromatography		
Endotoxin:	≤ 0.01 EU/µg by LAL analysis		
Storage & Stability:	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

