

## GLP1R

### Research Grade Anti-Human GLP1R (GMA105) mAb

<b>Catalog No.</b>	DHE40101A DHE40101B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	Glucagon-like peptide 1 receptor, GLP-1 receptor, GLP-1-R, GLP-1R, GMA105		
<b>Description:</b>	Research Grade Anti-Human GLP1R (GMA105) was affinity purified from cell culture supernatant.		
<b>UniProt ID (target):</b>	P43220		
<b>Specificity:</b>	Human Glucagon-like peptide 1 receptor		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Source:</b>	Mammalian cell line		
<b>Isotype:</b>	Human IgG2		
<b>Formulation:</b>	Sterile filtered 0.01M PBS buffer, pH 7.4		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Clone:</b>	GMA105		
<b>Purification:</b>	Protein A/G purified from cell culture supernatant		
<b>Applications:</b>	Functional studies		
<b>Endotoxin:</b>	≤ 0.01 EU/µg by LAL analysis		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. <b>Avoid freeze/thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

