

## KIR2DL1, 2DL2, 2DL3 Research Grade Lirilumab

<b>Catalog No.</b>	DHE43401A DHE43401B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	BMS-986015, IPH2102, CAS: 1000676-41-4		
<b>Description:</b>	Research Grade Lirilumab is a human monoclonal antibody which targets Killer cell immunoglobulin-like receptor 2DL1 & Killer cell immunoglobulin-like receptor 2DL2 & Killer cell immunoglobulin-like receptor 2DL3		
<b>UniProt ID (target):</b>	P43626, P43627, P43628		
<b>Specificity:</b>	Human Killer cell immunoglobulin-like receptor 2DL1, 2DL2, 2DL3		
<b>Source:</b>	XtenCHO cell line		
<b>Isotype:</b>	Human IgG4-kappa		
<b>INN:</b>	Lirilumab		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Formulation:</b>	0.01M PBS, pH 7.4		
<b>Purity:</b>	> 95% as determined by SDS-PAGE		
<b>Purification:</b>	Protein A/G affinity chromatography		
<b>Endotoxin:</b>	≤ 0.01 EU/µg by LAL analysis		
<b>Storage &amp; Stability:</b>	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. <b>Avoid freeze/thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.