

CCR5

Research Grade Leronlimab

Catalog No.	DHE85401A DHE85401B	Quantity:	100 µg 1.0 mg
Alternate Names:	HuPRO-140, PA-14, PRO-140, CAS: 674782-26-4		
Description:	Research Grade Leronlimab is a humanized monoclonal antibody which targets human C-C chemokine receptor type 5		
UniProt ID (target):	P51681		
Specificity:	Human C-C chemokine receptor type 5 (CD195)		
Source:	XtenCHO cell line		
Isotype:	Human IgG4 kappa		
INN:	Leronlimab		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	0.01M PBS, pH 7.4		
Purity:	> 95% as determined by SDS-PAGE		
Purification:	Protein A/G affinity chromatography		
Endotoxin:	≤ 0.01 EU/µg by LAL analysis		
Storage & Stability:	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

