

## FCGRT

### Research Grade Nipocalimab

|                                 |                                                                                                                                 |                  |                  |
|---------------------------------|---------------------------------------------------------------------------------------------------------------------------------|------------------|------------------|
| <b>Catalog No.</b>              | DHF14604A<br>DHF14604B                                                                                                          | <b>Quantity:</b> | 100 µg<br>1.0 mg |
| <b>Alternate Names:</b>         | M-281, M281, CAS: 2211985-36-1                                                                                                  |                  |                  |
| <b>Description:</b>             | Research Grade Nipocalimab is a human monoclonal antibody which targets human IgG receptor FcRn large subunit p51               |                  |                  |
| <b>UniProt ID (target):</b>     | P55899                                                                                                                          |                  |                  |
| <b>Target:</b>                  | Human IgG receptor FcRn large subunit p51                                                                                       |                  |                  |
| <b>Concentration:</b>           | 1.0 mg/ml, lot specific                                                                                                         |                  |                  |
| <b>Source:</b>                  | XtenCHO                                                                                                                         |                  |                  |
| <b>Isotype:</b>                 | Human IgG1 lambda                                                                                                               |                  |                  |
| <b>Formulation:</b>             | 0.01M PBS buffer, pH 7.4                                                                                                        |                  |                  |
| <b>Purity:</b>                  | > 95% by SDS-PAGE                                                                                                               |                  |                  |
| <b>Endotoxin Level:</b>         | ≤ 0.01 EU/µg by LAL analysis                                                                                                    |                  |                  |
| <b>Purification:</b>            | Protein A/G affinity chromatography                                                                                             |                  |                  |
| <b>INN:</b>                     | Nipocalimab                                                                                                                     |                  |                  |
| <b>Applications:</b>            | Functional studies                                                                                                              |                  |                  |
| <b>Storage &amp; Stability:</b> | Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year.<br><b>Avoid freeze/thaw cycles.</b> |                  |                  |

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



**Cell Sciences**<sup>®</sup>  
65 Parker Street  
Unit 11  
Newburyport, MA 01950

Toll Free: 888-769-1246  
Phone: 978-572-1070  
Fax: 978-992-0298

E-mail: [info@cellsciences.com](mailto:info@cellsciences.com)  
Website: [www.cellsciences.com](http://www.cellsciences.com)