

## CLDN18 & CD3E Research Grade Gresonitamab

<b>Catalog No.</b>	DHF19802A DHF19802B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	Bispecific, AMG 910, AMG-910, AMG910 CAS: 2413817-97-5		
<b>Description:</b>	Research Grade Gresonitamab is a bispecific human/chimeric/humanized monoclonal antibody that targets human Claudin-18 and T-cell surface glycoprotein CD3 epsilon chain		
<b>UniProt ID (target):</b>	P56856-2 & P07766		
<b>Specificity:</b>	Human Claudin-18 and T-cell surface glycoprotein CD3 epsilon chain		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Source:</b>	XtenCHO cell line		
<b>Isotype:</b>	Human (scFv-heavy-kappa)-(scFv-heavy-lambda)-scFc		
<b>Formulation:</b>	0.01M PBS buffer, pH 7.4		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Endotoxin Level:</b>	≤ 0.01 EU/µg by LAL analysis		
<b>Purification:</b>	Protein A/G affinity chromatography		
<b>Applications:</b>	Functional studies		
<b>INN:</b>	Gresonitamab		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. <b>Avoid freeze/thaw cycles.</b>		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**