

CXCR4

Research Grade Anti-Human CD184/CXCR4 (PF-06747143)

Catalog No.	DHF39909A DHF39909B	Quantity:	100 µg 1.0 mg
Alternate Names:	C-X-C chemokine receptor type 4, CXCR4, CD184, SDF-1 receptor, LAP-3, LPS-associated protein 3, PF-06747143, PF 06747143, PF06747143		
Description:	Humanized monoclonal antibody PF-06747143 targeting C-X-C chemokine receptor type 4, CXCR4. Antagonist antibody potently inhibits CXCR4-mediated signaling pathways.		
UniProt ID:	P61073		
Specificity:	Human CXCR4		
Concentration:	1.0 mg/ml, lot specific		
Source:	Mammalian cell line		
Isotype:	Human IgG		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com