

SIRPA

Research Grade Anti-Human CD172a/SIRPA (CC-95251)

Catalog No.	DHF64006A DHF64006B	Quantity:	100 µg 1.0 mg
Alternate Names:	Tyrosine-protein phosphatase non-receptor type substrate 1, SIRP-alpha, CD172a, CC-95251, CC 95251, CC95251		
Description:	Research Grade Anti-Human CD172a/SIRPA is a high affinity human monoclonal antibody which targets human Tyrosine-protein phosphatase non-receptor type substrate 1 (SIRP-alpha) and blocks the CD47- SIRPα interaction.		
UniProt ID:	P78324		
Specificity:	Human SIRP-alpha (CD172a)		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Human IgG		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A or G affinity chromatography		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com