

ADAM17

Research Grade Anti-Human CD156b/ADAM17 (D8P1C1)

Catalog No.	DHF68802A DHF68802B	Quantity:	100 µg 1.0 mg
Alternate Names:	Disintegrin and metalloproteinase domain-containing protein 17, ADAM 17, CD156b, Snake venom-like protease, TNF-alpha convertase, TNF-alpha-converting enzyme		
Description:	Research Grade Anti-Human CD156b/ADAM17 (D8P1C1) is a human monoclonal antibody which targets TNF-alpha-converting enzyme (ADAM 17, CD156b)		
UniProt ID (target):	P78536		
Specificity:	Human TNF-alpha-converting enzyme (ADAM 17, CD156b)		
Source:	XtenCHO cell line		
Isotype:	Human IgG		
Clone:	D8P1C1		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	0.01M PBS, pH 7.4.		
Purity:	> 95% as determined by SDS-PAGE		
Purification:	Protein A/G affinity chromatography		
Endotoxin:	≤ 0.01 EU/µg by LAL analysis		
Applications:	Functional studies		
Storage & Stability:	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

