

DLK1

Research Grade Anti-Human DLK1 (ADCT-701) mAb

Catalog No.	DHF71402A DHF71402B	Quantity:	100 µg 1.0 mg
Alternate Names:	Protein delta homolog 1, DLK1, pG2, FA1, ADCT-701, ADCT 701, ADCT701		
Description:	Research Grade Anti-Human DLK1 (ADCT-701) is a humanized monoclonal antibody which targets human protein delta homolog (DLK1)		
UniProt ID (target):	P80370		
Specificity:	Human protein delta homolog (DLK1)		
Source:	Mammalian cells		
Isotype:	Human IgG1 kappa		
Clone:	ADCT-701		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	0.01M PBS, pH 7.4.		
Purity:	> 95% as determined by SDS-PAGE		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Endotoxin:	≤ 0.001 EU/µg by LAL analysis		
Applications:	Functional studies		
Storage & Stability:	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

