

TNFRSF17

Research Grade Belantamab

Catalog No.	DHF92401A DHF92401B	Quantity:	100 µg 1.0 mg
Alternate Names:	BCMA, 2857916, GSK2857916, J6M0, J6M0-mcMMAF, CAS: 2061894-48-0		
Description:	Research Grade Belantamab is a humanized monoclonal antibody directed against Human Tumor necrosis factor receptor superfamily member 17 (CD269, BCMA)		
UniProt ID (target):	Q02223		
Target:	Human Tumor necrosis factor receptor superfamily member 17 (CD269, BCMA)		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Human IgG1 kappa		
Formulation:	Sterile-filtered 0.01 M PBS, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
INN:	Belantamab		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences®
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com