

## TNFRSF17 & CD3E Research Grade Pacanalotamab

<b>Catalog No.</b>	DHF92404A DHF92404B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	Bispecific, AMG-420, AMG420, BI-836909, CAS: 2251756-52-0		
<b>Description:</b>	Research Grade Pacanalotamab is a bispecific humanized monoclonal antibody directed against human Tumor necrosis factor receptor superfamily member 17 (BCMA, CD269) and T-cell surface glycoprotein CD3 epsilon chain (CD3e)		
<b>UniProt ID (target):</b>	Q02223 & P07766		
<b>Target:</b>	Human TNFRSF17 and CD3e		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Source:</b>	XtenCHO		
<b>Isotype:</b>	Human (scFv-heavy-kappa)-(scFv-heavy-lambda)-scFc		
<b>Formulation:</b>	30 mM histidine pH 5.8, 10% sucrose, 0.06% Tween 80		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Endotoxin Level:</b>	≤ 0.01 EU/µg by LAL analysis		
<b>Purification:</b>	Protein A/G affinity chromatography		
<b>INN:</b>	Pacanalotamab		
<b>Applications:</b>	Functional studies		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. <b>Avoid freeze/thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

