

FOLH1 & CD3E Research Grade Pasotuxizumab

Catalog No.	DHG01301A DHG01301B	Quantity:	100 µg 1.0 mg
Alternate Names:	Bispecific, BAY 2010112, BITE MT112, MT-112, CAS: 1442657-12-6		
Description:	Research Grade Pasotuxizumab is a bispecific chimeric humanized monoclonal which targets human Glutamate carboxypeptidase 2 and T-cell surface glycoprotein CD3 epsilon chain		
UniProt ID (target):	Q04609 & P07766		
Target:	Human Glutamate carboxypeptidase 2 and CD3E		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Chimeric Human (scFv-heavy-kappa)-(scFv-heavy-lambda)		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
INN:	Pasotuxizumab		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.