

FOLH1

Research Grade HuJ591-Gsmab

Catalog No.	DHG01302A DHG01302B	Quantity:	100 µg 1.0 mg
Description:	Research Grade HuJ591-Gsmab is a humanized monoclonal which targets human Glutamate carboxypeptidase 2 (Folate hydrolase 1)		
UniProt ID (target):	Q04609		
Target:	Human Glutamate carboxypeptidase 2		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Human IgG1		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
INN:	HuJ591-Gsmab		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

