

FOLH1

Research Grade Rosopatamab

Catalog No.	DHG01304A DHG01304B	Quantity:	100 µg 1.0 mg
Alternate Names:	J-591, MLN-591, huJ-591, HuJ591-DOTA, Rosopatamab tetraxetan, CAS: 2260767-49-3		
Description:	Research Grade Rosopatamab is a chimeric monoclonal antibody which targets human Glutamate carboxypeptidase 2 (Folate hydrolase 1)		
UniProt ID (target):	Q04609		
Target:	Human Glutamate carboxypeptidase 2		
Concentration:	1.0 mg/ml, lot specific		
Source:	Mammalian cells		
Isotype:	Chimeric IgG1 kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
INN:	Rosopatamab		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com