

## FOLH1 & CD3E Research Grade Voxalatamab

<b>Catalog No.</b>	DHG01307A DHG01307B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	Bispecific, JNJ-63898081, JNJ-8081 CAS: 2411871-58-2		
<b>Description:</b>	Research Grade Voxalatamab is a bispecific human/chimeric/humanized monoclonal antibody that targets human Glutamate carboxypeptidase 2 (GCP2) and T-cell surface glycoprotein CD3 epsilon chain		
<b>UniProt ID (target):</b>	Q04609 & P07766		
<b>Specificity:</b>	Human Glutamate carboxypeptidase 2 (GCP2) and T-cell surface glycoprotein CD3 epsilon chain		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Source:</b>	XtenCHO cell line		
<b>Isotype:</b>	Human Ig(G4 kappa_G4 lambda)		
<b>Formulation:</b>	0.01M PBS buffer, pH 7.4		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Endotoxin Level:</b>	≤ 0.01 EU/µg by LAL analysis		
<b>Purification:</b>	Protein A/G affinity chromatography		
<b>Applications:</b>	Functional studies		
<b>INN:</b>	Voxalatamab		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. <b>Avoid freeze/thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

