

CD47

Research Grade Magrolimab

Catalog No.	DHG17601A DHG17601B	Quantity:	100 µg 1.0 mg
Alternate Names:	CD47, Hu5F9-G4, CAS: 2169232-81-7		
Description:	Research Grade Magrolimab is a chimeric humanized monoclonal antibody which targets Human Leukocyte surface antigen CD47 (MER6)		
UniProt ID (target):	Q08722		
Specificity:	Human Leukocyte surface antigen CD47 (MER6)		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Chimeric/Human IgG4 kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A or G affinity chromatography		
Applications:	Functional studies		
INN:	Magrolimab		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

