

SLC39A6

Research Grade Ladiratuzumab

Catalog No.	DHG52901A DHG52901B	Quantity:	100 µg 1.0 mg
Alternate Names:	hLIV22, SGN-LIV1A, CAS: 1629760-28-6		
Description:	Research Grade Ladiratuzumab is a humanized monoclonal antibody which targets human Zinc transporter ZIP6		
UniProt ID (target):	Q13433		
Target:	Human Zinc transporter ZIP6		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Human IgG1 kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
INN:	Ladiratuzumab		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

