

## IL18

### Research Grade Anti-Human IL18/IL1F4 (GSK 1070806)

<b>Catalog No.</b>	DHG70704A DHG70704B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	GSK-1070806, GSK 1070806, GSK1070806		
<b>Description:</b>	Research Grade Anti-Human IL18/IL1F4 (GSK 1070806) is a humanized IgG1 monoclonal antibody that binds to human IL-18 with high affinity and neutralizes its function.		
<b>UniProt ID (target):</b>	Q14116		
<b>Specificity:</b>	Human Interleukin-18 (IL-18) / IFN-gamma-inducing factor (IL1F4)		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Source:</b>	XtenCHO		
<b>Isotype:</b>	Human IgG1		
<b>Formulation:</b>	0.01M PBS buffer, pH 7.4		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Endotoxin Level:</b>	≤ 0.01 EU/µg by LAL analysis		
<b>Purification:</b>	Protein A/G affinity chromatography from cell culture supernatant		
<b>Applications:</b>	Functional studies		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. <b>Avoid freeze/thaw cycles.</b>		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**