

IL11RA

Research Grade Anti-Human IL11RA (Hu16E12) mAb

Catalog No.	DHG84902A DHG84902B	Quantity:	100 µg 1.0 mg
Alternate Names:	Interleukin-11 receptor subunit alpha, IL-11 receptor subunit alpha, IL-11R subunit alpha, IL-11R-alpha; IL-11RA, IL-11R alpha, Hu16E12		
Description:	Human monoclonal antibody (clone Hu16E12) which targets human Interleukin-11RA		
UniProt ID (target):	Q14626		
Source:	Mammalian cells		
Specificity:	Human IL-11RA		
Isotype:	Human IgG4		
Clone:	Hu16E12		
Purification:	Protein A or G affinity chromatography		
Purity:	> 95% by SDS-PAGE		
Formulation:	Sterile-filtered 0.01M PBS, pH 7.4		
Endotoxin:	< 0.01 EU/µg		
Concentration:	1.0 mg/ml, lot specific		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for immediate use (1-2 weeks) or aliquot the stock solution and store at -20°C to -80°C up to 1 year. Avoid repeated freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

