

## PDCD1

### Research Grade Retifanlimab

<b>Catalog No.</b>	DHH02218A DHH02218B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	INCMGA-00012, MGA-012, CAS: 2079108-44-2		
<b>Description:</b>	Research Grade Retifanlimab is a chimeric humanized monoclonal antibody which targets human Programmed cell death protein 1 (PD-1, CD279)		
<b>UniProt ID (target):</b>	Q15116		
<b>Target:</b>	Human Programmed cell death protein 1 (PD-1, CD279)		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Source:</b>	XtenCHO		
<b>Isotype:</b>	Chimeric human IgG4 kappa		
<b>Formulation:</b>	0.01M PBS buffer, pH 7.4		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Endotoxin Level:</b>	≤ 0.01 EU/µg by LAL analysis		
<b>Purification:</b>	Protein A/G affinity chromatography		
<b>INN:</b>	Retifanlimab		
<b>Applications:</b>	Functional studies		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. <b>Avoid freeze/thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

