

PDCD1

Research Grade Nofazinlimab

Catalog No.	DHH02230A DHH02230B	Quantity:	100 µg 1.0 mg
Alternate Names:	CS 1003, CS-1003, CS1003 CAS: 2377845-98-0		
Description:	Research Grade Nofazinlimab is a humanized monoclonal antibody that targets human Programmed cell death protein 1 (CD279)		
UniProt ID (target):	Q15116		
Specificity:	Human Programmed cell death protein 1 (CD279)		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO cell line		
Isotype:	Human IgG4 kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
Applications:	Functional studies		
INN:	Nofazinlimab		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

