

PDCD1

Research Grade Zimberelimab

Catalog No.	DHH02235A DHH02235B	Quantity:	100 µg 1.0 mg
Alternate Names:	CAS: 2259860-24-5		
Description:	Research Grade Zimberelimab is a human monoclonal antibody which targets human programmed cell death 1 protein (PD-1) receptor.		
UniProt ID:	Q15116		
Specificity:	Human PD-1 (CD279)		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Human IgG4 lambda		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
Applications:	Functional studies		
INN:	Zimberelimab		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

