

IL17A

Research Grade Ixekizumab

Catalog No.	DHH28801A DHH28801B	Quantity:	100 µg 1.0 mg
Alternate Names:	LY2439821, CAS: 1143503-69-8		
Description:	Research Grade Ixekizumab is a humanized monoclonal antibody directed against Human Interleukin-17A		
UniProt ID (target):	Q16552		
Target:	Human Interleukin-17A		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Human IgG4 kappa		
Formulation:	Sterile-filtered 0.01 M PBS, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
INN:	Ixekizumab		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.