

IL17A & IL17F & ALB Research Grade Sonelokimab

Catalog No.	DHH28808A DHH28808B	Quantity:	100 µg 1.0 mg
Alternate Names:	Bispecific, Trivalent, M-1095, MSB-0010841, CAS: 1414386-05-2		
Description:	Research Grade Sonelokimab is a bispecific, trivalent, humanized monoclonal antibody which targets Human Interleukin-17A, Interleukin-17F and Albumin		
UniProt ID (target):	Q16552 & Q96PD4 & P02768		
Specificity:	Human IL-17A, IL-17F and albumin		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Human VH-VH'-VH		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A or G affinity chromatography		
Applications:	Functional studies		
INN:	Sonelokimab		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.