

## CD276 & CD3E

### Research Grade Obrindatamab

<b>Catalog No.</b>	DHJ04004A DHJ04004B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	Bispecific, MGD009, RES281 M1.1, hBRCA84D-hXR32-MP3 M1.1, CAS: 2069959-72-2		
<b>Description:</b>	Research Grade Obrindatamab is a bispecific humanized monoclonal antibody directed against human CD276 and CD3e		
<b>UniProt ID (target):</b>	Q5ZPR3 & P07766		
<b>Target:</b>	Bispecific human CD276 and CD3e		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Source:</b>	XtenCHO		
<b>Isotype:</b>	Human scFv with Domain Crossover,Kappa;Lambda		
<b>Formulation:</b>	Sterile-filtered 0.01 M PBS, pH 7.4		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Endotoxin Level:</b>	≤ 0.01 EU/µg by LAL analysis		
<b>Purification:</b>	Protein A/G affinity chromatography		
<b>INN:</b>	Obrindatamab		
<b>Applications:</b>	Functional studies		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. <b>Avoid freeze/thaw cycles.</b>		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**

