

CSPG4

Research Grade Tecnemab-K-1

Catalog No.	DHJ08301A DHJ08301B	Quantity:	100 µg 1.0 mg
Alternate Names:	225.28S, Tecnemab-K-1		
Description:	Research Grade Tecnemab-K-1 is a murine monoclonal antibody directed against human Chondroitin sulfate proteoglycan 4		
UniProt ID (target):	Q6UVK1		
Target:	Human Chondroitin sulfate proteoglycan 4		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Mouse IgG1 kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
INN:	Tecnemab-K-1		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.