

BTLA

Research Grade Anti-Human CD272/BTLA (ANB032)

Catalog No.	DHJ14403A DHJ14403B	Quantity:	100 µg 1.0 mg
Alternate Names:	ANB032		
Description:	ANB032 is a human monoclonal antibody directed against Human B- and T-lymphocyte attenuator (CD272), produced in XtenCHO cells.		
UniProt ID (target):	Q7Z6A9		
Specificity:	Human B- and T-lymphocyte attenuator (CD272)		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Human IgG4 kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

