

PCSK9

Research Grade Recaticimab

Catalog No.	DHJ24010A DHJ24010B	Quantity:	100 µg 1.0 mg
Alternate Names:	SHR-1209, CAS: 2361290-85-7		
Description:	Research Grade Recaticimab is a humanized monoclonal antibody directed against human Proprotein convertase subtilisin/kexin type 9 (PC9, NARC-1)		
UniProt ID (target):	Q8NBP7		
Target:	Human Proprotein convertase subtilisin/kexin type 9 (PC9, NARC-1)		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Human IgG1 kappa		
Formulation:	Sterile-filtered 0.01 M PBS, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
INN:	Recaticimab		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

