

CD200R1

Research Grade Ucenprubart

Catalog No.	DHJ27904A DHJ27904B	Quantity:	100 µg 1.0 mg
Alternate Names:	Cell surface glycoprotein CD200 receptor 1, CD200R1, CD200 cell surface glycoprotein receptor, Cell surface glycoprotein OX2 receptor 1, CAS: 2414559-57-0		
Description:	Research Grade Ucenprubart is a humanized monoclonal antibody which targets human cell surface glycoprotein CD200 receptor 1 (CD200R1)		
UniProt ID (target):	Q8TD46		
Specificity:	Human cell surface glycoprotein CD200 receptor 1 (CD200R1)		
Concentration:	1.0 mg/ml, lot specific		
Source:	Mammalian cell line		
Isotype:	Humanized IgG4 kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A or G affinity chromatography		
Applications:	Functional studies		
INN:	Ucenprubart		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

