

HAVCR2

Research Grade Surzebiclimab

Catalog No.	DHJ28403A DHJ28403B	Quantity:	100 µg 1.0 mg
Alternate Names:	BGB-A425, CAS: 2342597-90-2		
Description:	Research Grade Surzebiclimab is a humanized monoclonal antibody directed against human Hepatitis A virus cellular receptor 2 (CD366, TIMD-3)		
UniProt ID (target):	Q8TDQ0		
Specificity:	Human Hepatitis A virus cellular receptor 2 (CD366, TIMD-3)		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO cell line		
Isotype:	Humanized IgG1 kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
Applications:	Functional studies		
INN:	Surzebiclimab		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

