

LRRC15

Research Grade Samrotamab

Catalog No.	DHJ28801A DHJ28801B	Quantity:	100 µg 1.0 mg
Alternate Names:	ABBV-085, PR-1498487-MMAE, PR-1498487, Samrotamab vedotin, CAS:2052591-32-7		
Description:	Research Grade Samrotamab is a chimeric humanized monoclonal antibody that targets human Leucine-rich repeat-containing protein 15		
UniProt ID (target):	Q8TF66		
Target:	Human Leucine-rich repeat-containing protein 15		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Chimeric human IgG1 kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
INN:	Samrotamab		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

