

MUC16

Research Grade Igovomab

Catalog No.	DHJ31002A DHJ31002B	Quantity:	100 µg 1.0 mg
Alternate Names:	Indimacis 125, OC125 chelator DTPA, CAS: 171656-50-1		
Description:	Research Grade Igovomab is a murine monoclonal antibody which targets human tumor-associated antigen, CA-125.		
UniProt ID (target):	Q8WXI7		
Target:	Human Mucin-16 (CA-125)		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Mouse F(ab') ₂ -G1-nd		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
INN:	Igovomab		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.