

HTRA1

Research Grade Galegenimab

Catalog No.	DHJ32901A DHJ32901B	Quantity:	100 µg 1.0 mg
Alternate Names:	FHTR2163, RG-6147, RG6417, RO-7171009, RO7171009 CAS: 2403683-24-7		
Description:	Research Grade Galegenimab is a humanized monoclonal antibody that targets human Serine protease HTRA1		
UniProt ID (target):	Q92743		
Specificity:	Human Serine protease HTRA1		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO cell line		
Isotype:	Human IgG2 kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
Applications:	Functional studies		
INN:	Galegenimab		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com