

## TSLP

### Research Grade Tezepelumab

<b>Catalog No.</b>	DHJ36101A DHJ36101B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	AMG 157, AMG-157, MEDI-9929, MEDI9929, anti-TSLP-Receptor (TSLP-R, CrI2), CAS: 1572943-04-4		
<b>Description:</b>	Research Grade Tezepelumab is a human monoclonal antibody which targets Thymic stromal lymphopoietin		
<b>UniProt ID (target):</b>	Q969D9		
<b>Specificity:</b>	Human Thymic stromal lymphopoietin		
<b>Source:</b>	XtenCHO cell line		
<b>Isotype:</b>	Human IgG2-lambda		
<b>INN:</b>	Tezepelumab		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Formulation:</b>	0.01M PBS, pH 7.4.		
<b>Purity:</b>	> 95% as determined by SDS-PAGE		
<b>Purification:</b>	Protein A/G affinity chromatography		
<b>Endotoxin:</b>	≤ 0.01 EU/µg by LAL analysis		
<b>Storage &amp; Stability:</b>	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. <b>Avoid freeze/thaw cycles.</b>		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**

