

NECTIN4

Research Grade Enfortumab

Catalog No.	DHJ41001A DHJ41001B	Quantity:	100 µg 1.0 mg
Alternate Names:	AGS-22CE, AGS-22M, AGS-22M6E, unconjugated : AGS-22C3 or AGSM6, CAS: 1346452-25-2		
Description:	Research Grade Enfortumab is a human monoclonal antibody directed against Human Nectin-4		
UniProt ID (target):	Q96NY8		
Target:	Human Nectin-4		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Human IgG1 kappa		
Formulation:	Sterile-filtered 0.01 M PBS, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
INN:	Enfortumab		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

