

IL25

Research Grade Anti-Human IL25/IL17E (ABM125)

Catalog No.	DHJ55502A DHJ55502B	Quantity:	100 µg 1.0 mg
Alternate Names:	Interleukin-25, Interleukin-17E, IL-17E		
Description:	Research Grade Anti-Human Interleukin 25/IL17E (ABM125) is an affinity purified human monoclonal antibody.		
UniProt ID (target):	Q9H293		
Target:	Human Interleukin-25 (IL-17E)		
Concentration:	1.0 mg/ml, lot specific		
Source:	Mammalian cells		
Isotype:	Human IgG		
Formulation:	Sterile filtered 0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Clone:	ABM125		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

