

IL21

Research Grade Anti-Human IL21 (NNC-114-0005)

Catalog No.	DHJ59402A DHJ59402B	Quantity:	100 µg 1.0 mg
Alternate Names:	NNC0114-0005		
Description:	Research Grade Anti-Human IL21 (NNC-114-0005) is a human monoclonal antibody which targets human Interleukin-21		
UniProt ID (target):	Q9HBE4		
Specificity:	Human Interleukin-21		
Source:	XtenCHO cell line		
Isotype:	Human IgG		
Clone:	NNC-114-0005		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	0.01M PBS, pH 7.4.		
Purity:	> 95% as determined by SDS-PAGE		
Purification:	Protein A/G affinity chromatography		
Endotoxin:	≤ 0.01 EU/µg by LAL analysis		
Storage & Stability:	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.