

IL23A

Research Grade Tildrakizumab

Catalog No.	DHJ63102A DHJ63102B	Quantity:	100 µg 1.0 mg
Alternate Names:	MK-3222, SCH900222, CAS: 1326244-10-3		
Description:	Research Grade Tildrakizumab is a humanized monoclonal antibody which targets Human Interleukin-23 subunit alpha (IL-23p19)		
UniProt ID (target):	Q9NPF7		
Specificity:	Human Interleukin-23 subunit alpha (IL-23p19)		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Human IgG1 kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A or G affinity chromatography		
Applications:	Functional studies		
INN:	Tildrakizumab		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.