

## IL1RAP

### Research Grade Anti-Human IL1RAP (CAN10)

<b>Catalog No.</b>	DHJ63402A DHJ63402B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	Interleukin-1 receptor accessory protein, IL1RAP, Interleukin-1 receptor 3, IL-1R-3, IL-1R3		
<b>Description:</b>	Research Grade Anti-Human IL1RAP (CAN10) is a human monoclonal antibody which targets human Interleukin-1 receptor accessory protein		
<b>UniProt ID (target):</b>	Q9NPH3		
<b>Target:</b>	Human Interleukin-1 receptor accessory protein (IL1RAP)		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Source:</b>	Mammalian cells		
<b>Isotype:</b>	Human IgG1 kappa		
<b>Formulation:</b>	0.01M PBS buffer, pH 7.4		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Endotoxin Level:</b>	≤ 0.01 EU/µg by LAL analysis		
<b>Purification:</b>	Protein A/G affinity chromatography from cell culture supernatant		
<b>Clone:</b>	CAN10		
<b>Applications:</b>	Functional studies		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. <b>Avoid freeze/thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.