

IL20

Research Grade Fletikumab

Catalog No.	DHJ69101A DHJ69101B	Quantity:	100 µg 1.0 mg
Alternate Names:	15D2, NN-8226, NNC-0109-0012, CAS: 1357158-22-5		
Description:	Research Grade Fletikumab is a human monoclonal antibody which targets Interleukin -20		
UniProt ID (target):	Q9NYY1		
Specificity:	Human Interleukin-20		
Source:	XtenCHO cell line		
Isotype:	Human IgG4-kappa		
INN:	Fletikumab		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	0.01M PBS, pH 7.4.		
Purity:	> 95% as determined by SDS-PAGE		
Purification:	Protein A/G affinity chromatography		
Endotoxin:	≤ 0.01 EU/µg by LAL analysis		
Storage & Stability:	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.