

CD274

Research Grade Atezolizumab

| | | | |
|---------------------------------|---|------------------|------------------|
| Catalog No. | DHJ70101A DHJ70101B | Quantity: | 100 µg 1.0 mg |
| Alternate Names: | MPDL3280A, RG7446, CAS: 1380723-44-3 | | |
| Description: | Research Grade Atezolizumab is a humanized monoclonal antibody directed against Human Programmed cell death 1 ligand 1 | | |
| UniProt ID (target): | Q9NZQ7 | | |
| Target: | Human Programmed cell death 1 ligand 1 (CD274) | | |
| Concentration: | 1.0 mg/ml, lot specific | | |
| Source: | XtenCHO | | |
| Isotype: | Human IgG1 kappa | | |
| Formulation: | Sterile-filtered 0.01 M PBS, pH 7.4 | | |
| Purity: | > 95% by SDS-PAGE | | |
| Endotoxin Level: | ≤ 0.01 EU/µg by LAL analysis | | |
| Purification: | Protein A/G affinity chromatography | | |
| INN: | Atezolizumab | | |
| Applications: | Functional studies | | |
| Storage & Stability: | Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles. | | |

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences®
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com