

CD274

Research Grade Opucolimab

Catalog No.	DHJ70109A DHJ70109B	Quantity:	100 µg 1.0 mg
Alternate Names:	AHLX-20, HLX-20, CAS: 2251771-79-4		
Description:	Research Grade Opucolimab is a humanized monoclonal antibody which targets human Programmed cell death 1 ligand 1 (PD-L1, CD274)		
UniProt ID (target):	Q9NZQ7		
Specificity:	Human Programmed cell death 1 ligand 1 (CD274)		
Source:	XtenCHO cell line		
Isotype:	Human IgG1 lambda		
INN:	Opucolimab		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	0.01M PBS, pH 7.4		
Purity:	> 95% as determined by SDS-PAGE		
Purification:	Protein A/G affinity chromatography		
Endotoxin:	≤ 0.01 EU/µg by LAL analysis		
Storage & Stability:	Recommended storage, undiluted at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences®
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com